

# Quick Reference Guide

## – FL320

### Description


The FL320 is a aluminum lift for bariatric use with a weight capacity of 320 kg. It has electrical raising and lowering functions, base-width adjustment and adjustable mast height in 3 positions. Armrests add reassurance and stability for the user when gait training, toileting etc. Equipped with 4" front wheels and 5" lockable rear wheels, in addition to an easy access emergency stop button, the FL320 is an ideal floor lift for many types of lifting situations.



#### Part no:

FL320 kit North America: 71967-NA  
 FL320 kit Australia: 71967-AU  
 FL320 kit Great Britain: 71967-GB  
 FL320 kit Europe: 71967-EU

Kit includes: Lift, hand control, charger and 4-point hanger bar.



### Using the lift

<p>Select and position sling around the patient according to sling instructions.</p> <p>Use hand control to position carry bar directly above the patient.</p>	<p>Attach the head sling loops to carry bar first, than leg straps to leg hooks.</p>	<p>Raise the patient (up on hand control) and as raises check to ensure straps are firmly in place. Raise the patient up until obstacles are clear, hold sling and position/move patient as needed.</p>	<p>Lower patient (down on hand control) until in desired position. Continue to lower until all straps are slack to remove .</p>



**HUMAN CARE**  
 Your Life. Your Way.

**IMPORTANT!** This quick reference guide does not replace the user manual which can be downloaded at [www.humancaregroup.com](http://www.humancaregroup.com)

#### AUSTRALIA

Unit 9, 271 Wells Road  
 Chelsea Heights, Victoria 3196  
 Phone:- +61 3 8773 1111  
 Fax:- +61 3 9773 4008  
 Email: [info.au@humancaregroup.com](mailto:info.au@humancaregroup.com)

#### CANADA

10-155 Colonnade Road  
 Ottawa, ON K2E 7K1  
 Phone:- 613.723.6734  
 Fax:- 613.723.1058  
 Email: [info.ca@humancaregroup.com](mailto:info.ca@humancaregroup.com)

#### SWEDEN

Årstaängsvägen 21C  
 117 43 Stockholm  
 Phone:- +46 8 665 35 00  
 Fax:- +46 8 665 35 10  
 Email: [info.se@humancaregroup.com](mailto:info.se@humancaregroup.com)

#### UNITED STATES

4210 S. Industrial Drive, Suite 160  
 Austin, TX 78744  
 Phone:- 512.476.7199  
 Fax:- 512.476.7190  
 Email: [info.us@humancaregroup.com](mailto:info.us@humancaregroup.com)